


# PURELL SMARTLINK Activity Monitoring System Software Enhancements Quick Guide - July, 2021

New "Share Link" feature allows you to save your filtered report by creating a unique URL for future use or to share with other licensed users.

Reports now listed across the top. A few were discontinued due to low usage.

- EXPORT
- PDF/Print
- Excel
- Save Slide
- Raw Export



Expand / Hide Location Details allows you to see which floors, units and rooms you have included in the location filters for the report.



Edit Filters now includes all the filtering options in one area. All options are still available for timeframe. Shifts will appear if your facility has shifts set up in the software. All locations can be selected / unselected as needed.

